



## NATURAL CONVECTION OVEN BONC-105

## NATURAL CONVECTION OVEN BONC-105

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory natural Convection Oven.

## BONC-105 NATURAL CONVECTION OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain

Natural convection with low noise

PID controller with over temperature alarm and Timer function ensures precise and reliable control

## SPECIFICATIONS

Model	BONC-105
Capacity	200 L
Temperature Range	RT+10-200°C
Ambient Temperature	+5-35°C
Temperature Stability	±1°C
Display Resolution	0.1°C
Internal Dimension	570Wx560Dx640H mm
Exterior Dimension	710Wx780Dx945H mm
Shelves	3
Power Consumption	2200 W
Electrical Requirement	AC 220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)