



## NATURAL CONVECTION OVEN BONC-307

## NATURAL CONVECTION OVEN BONC-307

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Biochemical Studies, Hematological Studies, Bacterial Culturing, Research, Microbiological Determinations, Pharmaceutical Stability Assays, Food Processing Quality Control.

Also known as Laboratory Platelet Incubator.

## BONC-307 NATURAL CONVECTION OVEN



Microcomputer temperature controller

LCD Display

'Set temperature' and 'Actual temperature' display on double screen, clear and straight forward.

Air forced convection design, ensures the temperature uniform in the chamber.

Over temperature / Low temperature alarm, abnormal shutdown power outage alarm.

Sets in exchange stomatal, to ensure that the air is fresh in the cabinet.

UV Lamp in the cabinet.

Automatic defrost function.

Big observation window, working condition will be clear at a glance.

Stainless steel inner chamber material, easy to clean.

Heat preservation material, uses polyurethane foam.

Closed type compressor.

Less noise.

## SPECIFICATIONS

Model	BONC-307
Capacity	double doors, 20 floors for each sides, total 40 floors (80-120bags)
Temperature Range	15-35°C
Temperature Accuracy	±0.1
Timer Range	0-999
Variation (%)	±0.5
Internal Dimension	1090Wx390Dx950H mm
Exterior Dimension	1170Wx600Dx1750H mm
Tray Size (WxDxH)	980x300x810 mm
Weight	320 Kg
Power	1.35 KW



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)