



## NATURAL CONVECTION OVEN BONC-309

## NATURAL CONVECTION OVEN BONC-309

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry..

Also known as Vacuum Oven, Laboratory Vacuum Oven.

## BONC-309 NATURAL CONVECTION OVEN



- Stainless steel inner chamber without scroll rolled
- Silicone rubber door seal
- Strong sealing
- Nice design

## SPECIFICATIONS

Model	BONC-309
Capacity	25 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±2%
Vacuum Degree	< 133 Pa
Timer Range	0-999 min
Internal Dimension	300Wx300x280H mm
Shelves/Trays	2
Power	900 W
Power Supply	220 V, 50 Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)