



## NATURAL CONVECTION OVEN BONC-311

## NATURAL CONVECTION OVEN BONC-311

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry..

Also known as Vacuum Oven, Laboratory Vacuum Oven.

## BONC-311 NATURAL CONVECTION OVEN



- Stainless steel inner chamber without scroll rolled
- Silicone rubber door seal
- Strong sealing
- Nice design

## SPECIFICATIONS

Model	BONC-311
Capacity	50 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±2%
Vacuum Degree	< 133 Pa
Timer Range	0-999 min
Internal Dimension	415Wx350Dx400H mm
Shelves/Trays	2
Power	1400 W
Power Supply	220 V, 50 Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)