



HOLE-TYPE OR LAMP BOPL-201

HOLE-TYPE OR LAMP BOPL-201

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

BOPL-201 HOLE-TYPE OR LAMP



Combined type lighthouse

SPECIFICATIONS

| Model | BOPL-201 |
|----------------------------------------------------------------------|----------------------------|
| Illumination | 80,000LUX (at 1m distance) |
| Color temperature | 4500K±500K |
| Color rendering Index | 100≥Ra≥85 |
| Depth of light beam | >500 cm |
| The Value of total irradiance | <6 mW/m²lx |
| Total irradiance | <600 W/m² |
| Temperature rise on the head of the patient under surgical operation | ≤2°C |
| Rated Voltage of Bulb | 24V |
| Rated Power of bulb | 125 W |
| Input power | 150 VA |
| Power Supply | 220V±10% 50Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com