



HOLE-TYPE OR LAMP BOPL-203

HOLE-TYPE OR LAMP BOPL-203

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

BOPL-203 HOLE-TYPE OR LAMP



Color rendering Index (Ra) $100 \geq Ra \geq 85$

SPECIFICATIONS

Model	BOPL-203
Illumination	80,000LUX (at 1m distance)
Color temperature	4500K \pm 500K
Color rendering Index	$100 \geq Ra \geq 85$
Depth of light beam	>500 cm
The Value of total irradiance	<6 mW/m ² lx
Total irradiance	<600 W/m ²
Temperature rise on the head of the patient under surgical operation	$\leq 2^{\circ}\text{C}$
Rated Voltage of Bulb	24V
Rated Power of bulb	125 W
Input power	150 VA
Power Supply	220V \pm 10% 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com