



INTEGRATED OR LAMP BOPL-301

INTEGRATED OR LAMP BOPL-301

Shadowless operation lamp is suitable for the requirement of various operation conditions and it is ideal operating lighting equipment for modern operation rooms. Based on trichromatic theory, the light is closer to natural light, which is suitable for medical workers to distinguish the color of vascular tissue, nerve ending and so on. Adapts US blue cold light coating technology, the optical filter could effectively filter ultraviolet rays from the spectrum.

Used in Medical Urology, ENT, Gynecology examination, Surgeryunit, Surgery, Operation theatre.

Also known as Single-hole operaton lamp, Examination lamp, Surgical lamp, Operation light.

BOPL-301 INTEGRATED OR LAMP



SPECIFICATIONS

Model	BOPL-301
Illumination	60,000-120,000LX (at 1m distance)
Color temperature	4500K±500K
Color rendering Index	100≥Ra≥85
Diameter of lamp holder	500 mm
Shift of major and stand by lamps	automatic
Depth of light beam	≥700 cm
The Value of total irradiance	<6 mW/m ² lx
Total irradiance	<600 W/m ²
Brightness adjustment	CPU light- dimmer
Maximum Height for installation	2800-3200 mm
Light field diameter	10-28 cm at 1m
Bulb Power	150W 24V Bulb power
Input power	200+10% VA
Power Supply	220V±10% 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com