



INTEGRATED OR LAMP BOPL-302

INTEGRATED OR LAMP BOPL-302

Shadowless operation lamp is suitable for the requirement of various operation conditions and it is ideal operating lighting equipment for modern operation rooms. Based on trichromatic theory, the light is closer to natural light, which is suitable for medical workers to distinguish the color of vascular tissue, nerve ending and so on. Adapts US blue cold light coating technology, the optical filter could effectively filter ultraviolet rays from the spectrum.

Used in Urology, ENT, Gynecology examination, Surgery.

Also known as Single-hole operation lamp, Examination lamp, Surgical lamp, Operation light.

BOPL-302 INTEGRATED OR LAMP



Imported LED cold light source is adopted as operation lighting. As a genuine cold light source

Streamline design

The lamp can be flexibly regulated 180 degrees

Battery working (optional)

SPECIFICATIONS

Model	BOPL-302
Illumination	80,000-160,000LX (at 1m distance)
Color temperature	4500K±500K
Temperature rise	≤12℃
Temperature rise on the head of the patient	≤2℃
Diameter of lamp holder	500 mm
Shift of major and stand by lamps	automatic
Depth of light beam	≥600mm
Facula diameter	φ100~φ280 mm
Brightness adjustment	8-stage continuous
Bulb Voltage	AC24V Bulb power
Input Power	200VA
Power Supply	AC220V±10% 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com