



## LED OR SHADOWLESS LAMP BOPL-404

## LED OR SHADOWLESS LAMP BOPL-404

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

## BOPL-404 LED OR SHADOWLESS LAMP



Stepless adjusting light source  
ONDAL arm, OSRAM bulb

## SPECIFICATIONS

Model	BOPL-404
Illumination	40000~160000 LUX/40000~160000 LUX Illumination at 1m on Endo mode- 10000 LUX
Color temperature	6700>TC>3000K
Color rendering Index	100≥Ra≥85 Red color rendering Index R9-≥97
Light field diameter	D10 180+40/180+40 mm
Depth of L1+L2 illumination	150 cm
Heat to light ratio	<6mW/m²lx
D50	110+35/110+35
Number of LED bulbs	84/84
Average bulb life (hours)	50000 hours
Ceiling height for installation	270-320 cm
Input power	150W/150W VA
Power Supply	100-240V



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)