



DOUBLE ARM LED OPERATING LAMP BOPL-409

DOUBLE ARM LED OPERATING LAMP BOPL-409

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

BOPL-409 DOUBLE ARM LED OPERATING LAMP



Streamlined, ultra-thin, enclosed head with reduced air resistance to minimise contamination

Digital controlled light intensity without grades of brightness and automatic memory for added flexibility

Optionally built in power supply systems avoiding complicated installation of wall mounted control boxes

Auxiliary light source included for laparoscopic operations

Rigid suspension arm with easy mobility and stable positioning

No infrared or ultraviolet radiation

Backlight compensation function

Endoscopy light source

Stepless adjusting light source

SPECIFICATIONS

Model	BOPL-409
Illumination	160,000 Lux / 140,000 Lux
Color temperature	3600K-5500K (±500K) variable
Spot diameter (D50-D10)	130-350 / 130-350 mm
Color rendering Index	≥97 / ≥97
Average bulb life (hours)	100,000 h / 100,000 h Power of bulb 1 W x N 128 / 1 W x N 64 Light focus system-Manual/ Manual
Power Supply	220 V 50 Hz/ 220 V 50 Hz
Input power	160 W / 80 W
Ceiling height for installation	2.8 m / 2.8 m



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com