



ORTHOPEDICS TRACTOR RACK BOPT-401

ORTHOPEDICS TRACTOR RACK BOPT-401

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.

Used in Surgery, Leg operation, Medical unit.

Also known as Operating Table.

BOPT-401 ORTHOPEDICS TRACTOR RACK



Adopt carbon fiber material, the prespective imaging is clear, light and strong

The plug in design is convient and quick to connect with table

Equipped with a traction frame accessory car, which is convient to move and save space

SPECIFICATIONS

Model	BOPT-401
Table Length	1420±50 mm
Table Width	765±50 mm
Traction device horizontal rotary	0-180°
Traction distance	(0-200)+10%
Traction Device height adjustment	>400 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com