



## VACUUM OVEN BOVA-7203

## VACUUM OVEN BOVA-7203

## BOVA-7203 VACUUM OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time  
A big dual layer tempered glass on the door provides good observation  
Door adjustment system with silicon door gasket ensures better vacuity

## SPECIFICATIONS

Model	BOVA-7203
Capacity	24 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	700W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	300Wx300Dx275H mm
Exterior Dimension	605Wx490Dx450H mm
Shelves	2 pcs
Chamber Material	304 Stainless steel



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)