



VACUUM OVEN BOVA-7209

VACUUM OVEN BOVA-7209

BOVA-7209 VACUUM OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time
A big dual layer tempered glass on the door provides good observation
Door adjustment system with silicon door gasket ensures better vacuity

SPECIFICATIONS

Model	BOVA-7209
Capacity	30 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	550 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	320Wx320Dx300H mm
Exterior Dimension	630Wx510Dx490H mm
Shelves	1 pcs
Chamber Material	316 Stainless steel



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com