



VACUUM OVEN BOVA-7303

VACUUM OVEN BOVA-7303

BOVA-7303 VACUUM OVEN



Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

Direct heating and control by shelves: Each shelf has a separate temperature sensor, can be independently displayed and controlled temperature in setting, monitoring and timing, more accurate and stable

SPECIFICATIONS

Model	BOVA-7303
Capacity	215 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	2100 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	560Wx600Dx640H mm
Exterior Dimension	720Wx8200Dx1750H mm
Shelves	3 pcs
Chamber Material	304 Stainless steel



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com