



VACUUM OVEN BOVA-7306

VACUUM OVEN BOVA-7306

BOVA-7306 VACUUM OVEN



- Microprocessor controller
- Dual layer tempered glass door for clear observation
- Minimum heating time 50% less than traditional vacuum oven
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning
- Standard with vacuum pump

SPECIFICATIONS

| | |
|------------------------|---------------------|
| Model | BOVA-7306 |
| Capacity | 90 L |
| Electrical requirement | AC 220 V 50 HZ |
| Controller | LED display |
| Power Consumption | 1350 W |
| Temperature range | RT+10 - 200°C |
| Display Resolution | 0.1 |
| Temperature Stability | ±1°C |
| Vacuum Degree | 133 Pa |
| Ambient Temperature | +5 - 40°C |
| Interior Dimension | 450Wx450Dx450H mm |
| Exterior Dimension | 610Wx590Dx1350H mm |
| Shelves | 4 pcs |
| Chamber Material | 304 Stainless steel |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com