



VACUUM OVEN BOVA-7405

VACUUM OVEN BOVA-7405

BOVA-7405 VACUUM OVEN



- Microprocessor controller
- Dual layer tempered glass door for clear observation
- Minimum heating time 50% less than traditional vacuum oven
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning
- Standard with vacuum pump
- Provides programmability for 7 periods and 9 steps for each period making for a total of 63 programmable step

SPECIFICATIONS

Model	BOVA-7405
Capacity	63 L
Electrical requirement	AC 220 V 50 HZ
Controller	LCD Programmability controller
Power Consumption	1800 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	400Wx400Dx400H mm
Exterior Dimension	600Wx570Dx1390H mm
Shelves	3 pcs
Chamber Material	304 Stainless steel



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com