



PIPETTE PUMP BPIC-201

PIPETTE PUMP BPIC-201

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in variety of sizes and shapes.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Laboratory Pipette Pump.

BPIC-201 PIPETTE PUMP



- Economical
- Highly gas tightness
- Environment- friendly materials
- Lightweight and dexterous for operation
- Durable
- Recalibration Mechanism

SPECIFICATIONS

Model	BPIC-201
Volume	2ml
Color	Bule + White
Unit	1000



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com