



## PIPETTE STAND BPIC-301

## PIPETTE STAND BPIC-301

Pipette stand is a holder manufactured to provide safe and convenient storage for pipettes. Modern design enables this stand to securely hold virtually every pipettes in the marketplace, making it truly 'universal'.

Used in Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification.

Also known as Pipette Holder, Laboratory Pipette Stand, Laboratory Pipette Holder.

## BPIC-301 PIPETTE STAND



Linear type

Robust design

The freely rotating stand provides safe storage for all micropipettes

The stands are convenient to hold up 6 pipettes both of single and multi-channel.

## SPECIFICATIONS

Model	BPIC-301
Description	For 5 pipettes



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)