



## ULTRA LOW RETENTION FILTER TIP BPIC-416

## ULTRA LOW RETENTION FILTER TIP BPIC-416

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

## BPIC-416 ULTRA LOW RETENTION FILTER TIP



Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-416
Volume	10 $\mu$ l
Package	1000 pcs/bag
Unit	10



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)