



**RACKED ULTRA LOW RETENTION TIP 96X5
BPIC-424**

RACKED ULTRA LOW RETENTION TIP 96X5 BPIC-424

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

BPIC-424 RACKED ULTRA LOW RETENTION TIP 96X5



- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-424
Volume	300 µl
Unit	10



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com