



PIPETTOR SOLUTION BASINS BPIC-501

PIPETTOR SOLUTION BASINS BPIC-501

Basins are designed for single and multi channel pipettors. Easy access to solution. It avoids spills and provides stability. They are available pre-sterile or nonsterile and constructed of either polyvinyl chloride or polystyrene. Commonly used in laboratories, research, pharmaceutical, etc.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Laboratory Pipettor Solution Basins.

BPIC-501 PIPETTOR SOLUTION BASINS

Permits easy access to solution

Easy Cleaning

SPECIFICATIONS

Model	BPIC-501
Volume	50 ml



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com