

SEROLOGICAL PIPETTE BPIC-605


SEROLOGICAL PIPETTE BPIC-605

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

BPIC-605 SEROLOGICAL PIPETTE



Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation

SPECIFICATIONS

Model	BPIC-605
Volume	5 ml
Precision	2%
Negative Graduation (ml)	-2.5
Graduation	10-Jan
Color	Blue
Sterile	Gamma
Package	Individual
Pcs/bag	25 pcs/bag
Pcs/ctn	800 pcs/carton
Unit	1



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com