

## SEROLOGICAL PIPETTE BPIC-608


## SEROLOGICAL PIPETTE BPIC-608

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

## BPIC-608 SEROLOGICAL PIPETTE



Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation

## SPECIFICATIONS

Model	BPIC-608
Volume	10 ml
Precision	2%
Negative Graduation (ml)	-3
Graduation	10-Jan
Color	Red
Sterile	Gamma
Package	Bulk pack
Pcs/bag	25 pcs/bag
Pcs/ctn	3200 pcs/carton
Unit	1



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)