



STACKABLE 10?L TIP BOX BPIC-706

STACKABLE 10?L TIP BOX BPIC-706

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible.

Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

BPIC-706 STACKABLE 10?L TIP BOX



Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-706
Volume	10 μ l
Unit	20



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com