



10?L TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-712

10?L TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-712

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

BPIC-712 10?L TIP BOX WITH ULTRA LOW RETENTION TIP



Made from medical grade Polypropylene
Slippy surface with excellent design
DNase-free
Rnase-free
Non-pyrogenic
Filter and non-filter for choice
Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-712
Volume	10 µl
Unit	10



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com