



**200?L TIP BOX WITH ULTRA LOW RETENTION TIP
BPIC-713**

200µL TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-713

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

BPIC-713 200µL TIP BOX WITH ULTRA LOW RETENTION TIP



Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-713
Volume	200 µl
Unit	10



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com