



10000?L TIP BOX WITH ULTRA LOW RETENTION
TIP BPIC-717

10000µL TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-717

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

BPIC-717 10000µL TIP BOX WITH ULTRA LOW RETENTION TIP



- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-717
Volume	10000 µl
Unit	10



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com