



**300?L TIP BOX WITH ULTRA LOW RETENTION TIP  
BPIC-720**

## 300µL TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-720

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## BPIC-720 300µL TIP BOX WITH ULTRA LOW RETENTION TIP



Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-720
Volume	300 µl
Unit	10



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)