



## MOTORIZED PIPETTE CONTROLLER BPIP-413

## MOTORIZED PIPETTE CONTROLLER BPIP-413

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in laboratory and for research purpose.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Pipetaid, Pipet Gun, Laboratory Pipette Controller.

## BPIP-413 MOTORIZED PIPETTE CONTROLLER



- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling

## SPECIFICATIONS

|        |            |
|--------|------------|
| Model  | BPIP-413   |
| Volume | 0.1-100 ml |
| Colour | Dark Green |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)