



## PORTABLE REFRIGERATOR BPOR-100

## PORTABLE REFRIGERATOR BPOR-100

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## BPOR-100 PORTABLE REFRIGERATOR



## SPECIFICATIONS

Model	BPOR-100
External Size	330x235x275 mm
Internal Size	240x145x180 mm
Inner Capacity	6 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer
Insulation Material	PU
Gross Weight	2.8 kg
Packing Size	380x260x290 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)