



PORTABLE REFRIGERATOR BPOR-101

PORTABLE REFRIGERATOR BPOR-101

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

BPOR-101 PORTABLE REFRIGERATOR



SPECIFICATIONS

Model	BPOR-101
External Size	335x225x205 mm
Internal Size	295x167x168 mm
Inner Capacity	8 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer
Insulation Material	EPS
Gross Weight	1.8 kg
Packing Size	380x260x290 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com