



## PORTABLE REFRIGERATOR BPOR-102

## PORTABLE REFRIGERATOR BPOR-102

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## BPOR-102 PORTABLE REFRIGERATOR



## SPECIFICATIONS

Model	BPOR-102
External Size	425x262x305 mm
Internal Size	325x170x225 mm
Inner Capacity	12 L
Product Configuration	1pc refrigerator, 1pc tray 1pc partition 6 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU
Gross Weight	4 kg
Packing Size	450x330x330 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)