



PORTABLE REFRIGERATOR BPOR-103

PORTABLE REFRIGERATOR BPOR-103

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

BPOR-103 PORTABLE REFRIGERATOR



SPECIFICATIONS

Model	BPOR-103
External Size	425x262x305 mm
Internal Size	363x203x244 mm
Inner Capacity	15 L
Product Configuration	1pc refrigerator, 1pc strap 1pc tray 6 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU
Gross Weight	4.5 kg
Packing Size	450x330x330 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com