



## PORTABLE REFRIGERATOR BPOR-106

## PORTABLE REFRIGERATOR BPOR-106

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## BPOR-106 PORTABLE REFRIGERATOR



## SPECIFICATIONS

Model	BPOR-106
External Size	640x395x395 mm
Internal Size	590x335x320 mm
Inner Capacity	55 L
Product Configuration	1pc refrigerator, 12 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU
Gross Weight	6.5 kg
Packing Size	650x410x410 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)