



## PORTABLE REFRIGERATOR BPOR-204

## PORTABLE REFRIGERATOR BPOR-204

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## BPOR-204 PORTABLE REFRIGERATOR



1. The storage container can keep content cooling up to 24 hours
2. Insulation material is PU CFC-free polyurethane
3. Made of environmental PP material, food grade, non-toxic
4. Storage container is durable for vaccine storage, blood products storage, dry ice storage, frozen specimen storage etc

## SPECIFICATIONS

Model	BPOR-204
External Size (WxDxH)	570x310x370 mm
Internal Size (WxDxH)	450x225x300 mm
Capacity	33 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	20 h
Material	
External	PP
Internal	PU
Ice Vessels	3(CMB-1200)
Max Loading	10 kg
Optional Accessory	LED display thermometer
Net Weight	4.5 kg
Gross Weight	9.5 kg
Package Size (WxDxH)	590x330x390 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)