



LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER BPTL-201

LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER BPTL-201

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-201 LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER



Automatic calibration procedure, graphic display. Multi-point temperature calibration, can be added by user. It has tube locking function, can avoid influencing accuracy when tube is moved for high viscosity.

SPECIFICATIONS

Model	BPTL-201
Lifting method	One-click, Automatic
Viscosity and Gelation curve	Automatic scan
Gelation index and Temperature analysis	Automatic
Temperature	Multi-point calibration



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com