



SHEAR STABILITY TESTER BPTL-202

SHEAR STABILITY TESTER BPTL-202

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-202 SHEAR STABILITY TESTER



Independent desktop controller, high integrated function button easy operation. Automatic stroke counter can cut off power automatically.

SPECIFICATIONS

Model	BPTL-202
Timing method	Automatic stroke timer
Temperature control range	30 - 35 °C
Rotation speed	925 ± 25 r/min
Power	1100 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com