



CLOSED CUP FLASH POINT TESTER BPTL-207

CLOSED CUP FLASH POINT TESTER BPTL-207

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-207 CLOSED CUP FLASH POINT TESTER



Automatic forced air cooling after test. Flame size is adjustable. Protect program and self-diagnosis function.

SPECIFICATIONS

Model	BPTL-207
Lifting system	Automatic
Opening lid and Ignition	Automatic
Flame size	Adjustable
Air cooling	Automatic forced air cooling after test.
Touch screen	Displays test temperature,number of oil sampler,test time etc
Function	Protect program and self-diagnosis
Corrective function	Date,Time,Temperature and Atmospheric pressure



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com