



RUST PREVENTION TESTER BPTL-208

RUST PREVENTION TESTER BPTL-208

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-208 RUST PREVENTION TESTER



Independent quiet stirrer, can test one or several samples. Each stirrer controls rotating speed, test time and automatically stop when test is finished.

SPECIFICATIONS

Model	BPTL-208
Working units	4
Stirrer	Independent quiet, controls rotating speed, test time and automatically stop
Liquid circulator	Sputtering multi-channel
Mixing blade	Sliding automatic
Safety protection device	Prevent over-low liquid level



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com