



## COPPER STRIP CORROSION TESTER BPTL-211

## COPPER STRIP CORROSION TESTER BPTL-211

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-211 COPPER STRIP CORROSION TESTER



Can test lubricating oil and grease with different cylinders. Eco-friendly metal bath heating, no need heat-conducting medium. There is prompt when test is finished.

## SPECIFICATIONS

Model	BPTL-211
Samples size	6, simultaneously
Test samples	Lubricating oil and grease with different cylinders
Safety	Safety protection device to prevent over-high liquid level
Metal bath	Eco-friendly heating, no need heat-conducting medium



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)