



## WATER SEPARABILITY TESTER BPTL-212

## WATER SEPARABILITY TESTER BPTL-212

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-212 WATER SEPARABILITY TESTER



Sliding automatic location device, no need to adjust blade. Sputtering scroll liquid circulator can guarantee temperature uniformity of each point very well.

## SPECIFICATIONS

Model	BPTL-212
Working units	9
Stirrer	Digital, rotating, adjustable speed
Stirring time and record separation time	Automatically controlled
Unit of separation time	Selective
Safety	Safety protection device to prevent over-low liquid level
Liquid circulator	Sputtering scroll



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)