



DENSITY TESTER BPTL-213

DENSITY TESTER BPTL-213

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-213 DENSITY TESTER



Low noise magnetic pump cross circulation structure. Can test density of antifreeze

SPECIFICATIONS

Model	BPTL-213
Working units	2
Test sample	Antifreeze
Temperature control range	Ambient to 20 °C ± 0.1 °C
Refrigeration	Shock absorption compressor, fluorine free
Bath	Low temperature with lighting
Magnetic pump	Low noise with cross circulation structure
Safety	Safety protection device an prevent over-low liquid level



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com