



## DENSITY TESTER BPTL-213

## DENSITY TESTER BPTL-213

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-213 DENSITY TESTER



Low noise magnetic pump cross circulation structure. Can test density of antifreeze

## SPECIFICATIONS

|                           |   |
|---------------------------|---|
| Model                     | BPTL-213  |
| Working units             | 2   |
| Test sample               | Antifreeze  |
| Temperature control range | Ambient to 20 °C ± 0.1 °C                                 |
| Refrigeration             | Shock absorption compressor, fluorine free                |
| Bath                      | Low temperature with lighting                             |
| Magnetic pump             | Low noise with cross circulation structure                |
| Safety                    | Safety protection device an prevent over-low liquid level |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)