



Product Image Coming Soon

COOLING PERFORMANCE TESTER BPTL-218

COOLING PERFORMANCE TESTER BPTL-218

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-218 COOLING PERFORMANCE TESTER



Product Image Coming Soon

Curves of multi test results can be superimposed on one picture to contrast and analyze. Display maximum cooling speed, temperature of maximum cooling speed, cooling speed at 300°C, time etc.

SPECIFICATIONS

Model	BPTL-218
Probe	Released and lifted automatically
Cooling speed	300 °C
Test medium	Oil, Water, Quench liquid, Saline-alkaline solution, Steel rolling coolant, Cutting fluid.
Results	Curves of multi test results can be superimposed on one picture to contrast and analyze.
Display	Maximum cooling speed, Temperature of maximum cooling speed



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com