



## POUR POINT TESTER BPTL-221

## POUR POINT TESTER BPTL-221

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-221 POUR POINT TESTER



Dry type over speed refrigeration technology without medium. Temperature of two troughs are controlled separately. Compact structure, can put on platform directly

## SPECIFICATIONS

Model	BPTL-221
Trough	Two, temperature controlled separately
Samples size	2 samples in each trough
Temperature Range	Reach -30°C in 30 minutes and the lowest temperature is -70°C ,
Refrigeration	Dry type over speed refrigeration technology without medium
Refrigerant	Fluorine free
Structure	Compact structure, can put on platform directly



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)