



# CORROSION-PREVENTIVE PROPERTIES TESTER BPTL-223

## CORROSION-PREVENTIVE PROPERTIES TESTER BPTL-223

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-223 CORROSION-PREVENTIVE PROPERTIES TESTER



Bearing seat made of engineering plastics can avoid corrosive influence caused by metal bearing seat. Automatically stop motor when test is finished. Can do static and dynamic tests.

## SPECIFICATIONS

Model	BPTL-223
Tests	Static and Dynamic
Motor	Automatic stop when test is finished
Operation	Observe bearing through reflector with zoom function
Time	Can edit working time and stop time
Bearing seat	Made of engineering plastics can avoid corrosive influence



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)