



CORROSION TESTER BPTL-224

CORROSION TESTER BPTL-224

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-224 CORROSION TESTER



Sputtering multi-channel liquid circulator can guarantee temperature uniformity very well. Condenser quick detachable device

SPECIFICATIONS

Model	BPTL-224
Working units	3, can be extended to 4-6 units
Circulator	Sputtering multi-channel liquid circulator can guarantee temperature uniformity
Condenser	Quick detachable device
Safety protection device	Can prevent over-low liquid level



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com