



## VAPOUR PRESSURE TESTER BPTL-230

## VAPOUR PRESSURE TESTER BPTL-230

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-230 VAPOUR PRESSURE TESTER

Equipped with quick connector, make sure the test can be performed quickly and easily. Automatic temperature control, alarming while reaching required test temperature.

## SPECIFICATIONS

Model	BPTL-230
Test units	Two test units allow independent test for two samples at the same time
Display	Big colorful touch screen is convenient to use
Test accuracy	360° horizontal flip and oscillation improves test accuracy
Temperature control	Automatic temperature control, alarming while reaching required test temperature
Connector	Equipped with quick connector, make sure the test can be performed quickly and easily



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)