



NANO SPECTROPHOTOMETER BSNA-203

NANO SPECTROPHOTOMETER BSNA-203

Biolab Nano spectrophotometer is a compact, micro volume UV Visible Spectrophotometer with enhanced sensitivity that can detect samples at 3ng/μl. It gives you fast measurement typically less than 5sec/sample and provides rapid and reproducible results without prior dilution. With robust instrumentation and ease of use, our product is most preferred for quick analysis in labs.

Used in Ideal choice for quantitation of protein, DNA, RNA and photometric measurements in analytical labs and research areas..

BSNA-203 NANO SPECTROPHOTOMETER



Wavelength range: 200nm-800nm

Wavelength accuracy: ± 1.0 nm

Easy user interface

Reduce sample loss and eliminate the need for dilution by using low volumes of samples in 0.5μl to 2.0μl for sample measurement. The appropriate path length can be selected automatically or by manual selection.

Rapid measurements. Read-time is typically takes less than 5sec per sample without washing cuvettes or dilutes the samples

Reliability and robust instrumentation. Easy to use interface

Press-to-read feature reduces the amount of time .Optics with no moving parts minimizes the incidence of optical misalignment

Calibration curves, kinetics or ratio measurements are displayed at the touch of a button at connected computer.

Equipped with a built in printer

SPECIFICATIONS

Model	BSNA-203
Minimum Sample Size	0.3-2.0 μl
Wavelength Range	200-850 nm
Wavelength Accuracy	± 1 nm
Wavelength Reproducibility	$< \pm 0.2$ nm
Spectral Bandwidth	≤ 5 nm
Spectral Resolution	2 nm (FWHM at Hg 546)
Absorbance Range	0.002-75 (150, 300 optional, 10mm equivalent absorbance)
Absorbance Accuracy	1% (0.75 Abs at 350nm)
Absorbance Precision	0.002 Abs
Path Length	1mm, 0.2mm, (0.1mm, 0.05mm optional)
Light Source	Xenon flash lamp
Life Time of Lamp	10^{9+} flashes upto 10 years
Detector	2048-element linear silicon CCD array
Detection Range	Nucleic acid up to: 0.4~3750ng/μL(7500,15000 3750ng/μL(7500,15000 optional, dsDNA) and Protein up to: 0.01mg/ml~100mg/ml(200,400 optional, BSA)
Start-up Time	< 5 s
Measurement Time	< 5 s
Dimension (W/D/H)	240x210x110 mm
Weight (Net/Gross)	1.92 kg
Power	12V 4A 10 W
Power Supply	90-250V 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com