



## DECOLORIZATION SHAKER BSOT-205

## DECOLORIZATION SHAKER BSOT-205

Our product is an accurately designed, microprocessor controlled with double Decker platform to save your valuable lab space. Three eccentric shaft balancing drive ensures shaking with uniform speed. Speed adjustment settings permit both gentle and vigorous shaking. Varieties of platforms are available for different glassware and vessels. Also known as Laboratory Orbital Shaker.

## BSOT-205 DECOLORIZATION SHAKER



Compact appearance, space saving, powerful function and strong compatibility.  
Suction cup type machine feet, super shockproof, high speed of range 0~300 rpm and smooth, no noise.  
large platform design, placing more samples, strong load capacity.

## SPECIFICATIONS

Model	BSOT-205
Speed Range	0~300 rpm
Radius of gyration	10 mm
Max Load Capacity	3 Kgs
Motor Parameters	DC motor stepless speed regulation
Platform Dimension	330x430 mm
Beep Alarm	No
Voltage	AC110/220 V,50/60 Hz
Power	25 W
Fuse	250V,1A,φ5x20
Overall Dimension	W.330xD.490xH.145 mm
Net Weight	9.5 kgs

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
2001410006	Platform A	flat platform with rubber mat
2001410007	Platform B	platform with clamps for flasks 250-300 ml*12
2001410008	Platform C	platform with clamps for flasks 100-150 ml*20
2001410009	Platform D	universal platform with springs



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)